

# HoribaLabViewSDK

Antidoc v2.0.5, Niko Naredi-Rainer

# Table of Contents

1. Project description .....	1
2. DQMH® modules .....	2
2.1. Preamble .....	2
2.2. Modules overview .....	3
2.3. DeviceManager.lvlib .....	4
3. Classes .....	11
3.1. Classes overview .....	11
3.2. Communicator.lvclass .....	12
3.3. GenericDevice.lvclass .....	14
3.4. Monochromator.lvclass .....	16
3.5. CCD.lvclass .....	22
3.6. SCD.lvclass .....	37
4. Custom errors .....	39
5. Legal Information .....	40
5.1. Document creation .....	40
5.2. Product used in the project .....	42

# Chapter 1. Project description

The HORIBA SDK is a powerful tool for anyone working with HORIBA scientific instruments. It simplifies the process of instrument control and data acquisition, allowing users to focus on their research and development tasks. Whether you are conducting experiments, performing quality control, or developing new applications, this SDK provides the necessary tools to streamline your workflow and enhance productivity.

# Chapter 2. DQMH® modules

This section describes DQMH® module responsibilities and relationships.

## 2.1. Preamble

A DQMH module is the main component of an architecture based on DQMH® framework. A DQMH module is used to implement a section of the application that has one responsibility.

DQMH® framework defines two different type of DQMH module.

### **Singleton:**

A Singleton DQMH module can have only one instance running at any given time.

### **Cloneable:**

A Cloneable DQMH module can have one or multiple instances running in parallel.

DQMH® framework defines two different ways to carry data throughout the application and with both other DQMH modules and non-DQMH based code.

### **Request events:**

A request is a code that fires an event requesting the DQMH module to do something. Multiple locations in the code can send events to the DQMH module.

Request events are many-to-one.

Requests are usually named using imperative tense.

### **Broadcast events:**

A broadcast is a code that fires an event broadcasting that the DQMH module did something. Multiple Event Structures can register to handle the Broadcast Events.

Broadcast Events are one-to-many.

Broadcasts are usually named using past tense or passive voice.

### **NOTE**


Refer to the DQMH® framework official [documentation](#) to find more details on how the framework works

The following section gives you details on the project architecture relying on this framework. It


gives you an overview of the modules' interaction and detailed information on each module.

Graphs used in this section have the following legend:

**Components:**




**Events:**



**NOTE** One arrow can represent one or more events between two components

**NOTE** Request and Request and wait for Reply are represented by only one arrow. If there is no Request and wait for Reply, Request representation is used. Otherwise Request and wait for Reply is used

**Start and Stop module callers:**




## 2.2. Modules overview

This project contains 1 singleton module and 0 cloneable module.

*Table 1. Modules list*

Singleton	Cloneable
<a href="#">DeviceManager.lvlib</a>	

This graph represents the links between all DQMH modules.








## 2.3. DeviceManager.lvlib

**Type:** Singleton

**Responsibility:** This DQMH module takes care of starting and stopping the ICL as well as monitoring any output from the ICL. It also handles any communication to and from the ICL via Websocket requests.


### 2.3.1. Event list

*Table 2. Events*

Name	Type	Connector pane	Description	S.	R.	I.
Start Module			<p>Launches the module Main VI. After calling this VI, you can optionally register for broadcast events from the module by wiring the broadcast events output of this VI to a <b>Register For Events</b> function.</p> <p>After the optional Register For Events function call, you should always call the <b>Synchronize Module Events.vi</b> for this module with the 'Wait for Event Sync?' output of this VI to the corresponding input of the Synchronize Module Events.vi.</p> <p>To see an example of the proper wiring pattern, see the "Start Module: Value Change" event frame in the API Tester VI for this module.</p>			
Stop Module			<p>Send the Stop request to the Module's Main.vi.</p> <p>If <b>Wait for Module to Stop?</b> is TRUE, this VI will wait until the module main VI stops, and will timeout at the <b>Timeout to Wait for Stop</b> value. This value defaults to "-1", which means the VI will not timeout, and will always wait until the module main VI stops before completing execution.</p> <p>Note: The <b>Timeout to Wait for Stop</b> value is ignored if 'Wait for Module to Stop?' is set to FALSE.</p>			
Show Panel	→		Send the Show Panel request to the Module's Main.vi.			
Hide Panel	→		Send the Hide Panel request to the Module's Main.vi.			
Get Module Execution Status	→		Fire the Get Module Execution Status request.			

Name	Type	Connector pane	Description	S.	R.	I.
Show Diagram			This VI tells the Module to show its block diagram to facilitate troubleshooting (add probes, breakpoints, highlight execution, etc).			
SendJSONandAskForReply			This request asks in a synchronous manner to send a JSON string to the ICL, waits the specified time in [ms] and queries the websocket for a reply.			
StartICL			This command starts the ICL.exe and its monitoring			
OpenWebSocketCommunication			This event opens the websocket communication from the DeviceManager to the ICL.exe			
ICLshutdown			This request sends the command to shutdown the ICL.exe via websocket communication.			
DiscoverDevices			Requests from the ICL to discover monochromators, cameras and single channel detectors.			
DevicesList			This event calls mono_list, ccd_list, and scd_list.			
Module Did Init			Send the Module Did Init event to any VI registered to listen to this module's broadcast events.			
Status Updated			Send the Status Updated event to any VI registered to listen to events from the owning module.			
Error Reported			Send the Error Reported event to any VI registered to listen to events from the owning module.			
Module Did Stop			Send the Module Did Stop event to any VI registered to listen to this module's broadcast events.			
Update Module Execution Status			Broadcast event to specify whether or not the module is running.			



Name	Type	Connector pane	Description	S.	R.	I.
ICLstartNotification			This private event is used to tell the ICLcommunication loop that the ICL is running and a communication via websocket can be established			

**Type:**  → Request |  → Request and Wait for Reply |  → Broadcast

**Scope:**  → Protected |  → Community

**Reentrancy:**  → Preallocated reentrancy |  → Shared reentrancy

**Inlining:**  → Inlined

### 2.3.2. Module relationship




Table 3. Requests callers

Request Name	Callers
DevicesList	DeviceManager.lvlib:Test DeviceManager API.vi test_DevicesList.vi
DiscoverDevices	DeviceManager.lvlib:OpenConnectionWithDevice.vi DeviceManager.lvlib:Test DeviceManager API.vi test_DevicesDiscover.vi

Request Name	Callers
Get Module Execution Status	DeviceManager.lvlib:Obtain Broadcast Events for Registration.vi DeviceManager.lvlib:Start Module.vi
Hide Panel	DeviceManager.lvlib:Test DeviceManager API.vi
ICLshutdown	DeviceManager.lvlib:CloseConnectionWithDevice.vi DeviceManager.lvlib:Test DeviceManager API.vi
ICLstartNotification	DeviceManager.lvlib:Main.vi
OpenWebSocketCommunication	DeviceManager.lvlib:OpenConnectionWithDevice.vi DeviceManager.lvlib:Test DeviceManager API.vi test_DeviceManager_OpenWebSocketCommunication.vi
SendJSONandAskForReply	DeviceManager.lvlib:Test DeviceManager API.vi GenericDevice.lvclass:Send Receive Parse.vi
Show Diagram	DeviceManager.lvlib:Test DeviceManager API.vi
Show Panel	DeviceManager.lvlib:Test DeviceManager API.vi
StartICL	DeviceManager.lvlib:OpenConnectionWithDevice.vi DeviceManager.lvlib:Test DeviceManager API.vi test_DeviceManager_StartICL.vi

Table 4. Broadcasts Listeners

Broadcast Name	Listeners
Error Reported	DeviceManager.lvlib:Test DeviceManager API.vi
Module Did Init	DeviceManager.lvlib:Test DeviceManager API.vi
Module Did Stop	DeviceManager.lvlib:Test DeviceManager API.vi
Status Updated	DeviceManager.lvlib:Test DeviceManager API.vi
Update Module Execution Status	DeviceManager.lvlib:Test DeviceManager API.vi

Table 5. Used requests

Module	Requests
DeviceManager.lvlib	ICLstartNotification.vi Stop Module.vi

Table 6. Registered broadcast

Module	Broadcasts
DeviceManager.lvlib	Error Reported.vi Module Did Init.vi Module Did Stop.vi Status Updated.vi Update Module Execution Status.vi

### 2.3.3. Module Start/Stop calls




Table 7. Start and Stop module callers

Function	Callers
Start Module	DeviceManager.lvlib:OpenConnectionWithDevice.vi DeviceManager.lvlib:Test DeviceManager API.vi test_DeviceManager_Setup.vi
Stop Module	DeviceManager.lvlib:Handle Exit.vi DeviceManager.lvlib:CloseConnectionWithDevice.vi DeviceManager.lvlib:Test DeviceManager API.vi test_DeviceManager_Teardown.vi

### 2.3.4. Module custom errors

**TIP** Custom errors are added to the module via vi named **\*--error.vi**.

Module DeviceManager.lvlib use the following custom errors:

Table 8. Custom errors

Name	Code	Description
Module Not Running	0	
Module Not Stopped	0	
Module Not Synced	0	

Name	Code	Description
Request and Wait for Reply Timeout	0	

# Chapter 3. Classes


This section describes the classes contained in the project.

## 3.1. Classes overview

This project contains 5 classes and 0 interface.

*Table 9. Classes list*

Classes	Interfaces
<a href="#">Communicator.lvclass</a>	
<a href="#">GenericDevice.lvclass</a>	
<a href="#">Monochromator.lvclass</a>	
<a href="#">CCD.lvclass</a>	
<a href="#">SCD.lvclass</a>	



## 3.2. Communicator.lvclass

**Responsibility:** This class handles the communication between LV and the ICL via websocket requests.







Table 10. Functions (non private scope only)

Name	Connector pane	Description	S.	R.	I.
OpenConnect ion		Opens a connection to the ICL via websocket.			
CloseConnect ion		Closes the websocket connection.			
SendAndAsk ForReply		Wrapper around send and receive for the websocket communication.			
JSON Parse Command		Parses the original command.			
JSON Parse Errors		Parses the returned error into an LV error.			
JSON Parse ID		Parses the message ID to follow communication.			
JSON Parse Info		Wrapper around all other VIs that parse an answer from the ICL.			
JSON Parse Results		Parses the results from the ICL.			
Overwrite ID		Overwrites the message ID if a custom ID is to be used for messaging to the ICL.			

Name	Connector pane	Description	S.	R.	I.
ReceiveJson		Receives the reply from the synchronous communication to the ICL.			
SendJson		Sends the request for a synchronous communication to the ICL.			

Scope:  → Protected |  → Community

Reentrancy:  → Preallocated reentrancy |  → Shared reentrancy

Inlining:  → Inlined

### 3.3. GenericDevice.lvclass

**Responsibility:** This parent class for all devices provides base functionality like discovering, opening and closing devices on the ICL.

**Version:** 1.0.0.7







Table 11. Functions (non private scope only)

Name	Connector pane	Description	S.	R.	I.
DeviceClose		Closes communications with the CCD indicated by the index.			
DeviceOpen		This command initializes the CCD and gets it's the CCD configuration from the device. The device is also connected to the API. Since a CCD hardware initialization occurs, all CCD parameters, including any previously set parameters, will be reset to their default values.			
DeviceIsOpen		Returns true if selected CCD is open.			



Name	Connector pane	Description	S.	R.	I.
Send Receive Parse		Wrapper around the send and receive calls.			
ReadDeviceJSONList		Parses the returned devies into a string array.			
DecodeDeviceJSONList		Decodes the answer from the ICL into a device.			
ReadDeviceTypeAsString		Return the device class name as string.			
Read isOpen		Accessor VI for this class property.			
Write isOpen		Accessor VI for this class property.			
Read DeviceTerm		Accessor VI for this class property.			
Write DeviceTerm		Accessor VI for this class property.			
Read DeviceType		Accessor VI for this class property.			
Write DeviceType		Accessor VI for this class property.			
Read Index		Accessor VI for this class property.			
Write Index		Accessor VI for this class property.			
Read productId		Accessor to the product ID.s			
Write productId		Accessor to the product ID.			
Read SerialNumber		Accessor VI for this class property.			
Write SerialNumber		Accessor VI for this class property.			

Scope: → Protected | → Community

Reentrancy: → Preallocated reentrancy | → Shared reentrancy

Inlining: → Inlined

## 3.4. Monochromator.lvclass

**Responsibility:** This class contains all functionality needed for Horiba's Monochromators.

**Version:** 1.0.0.1




Table 12. Functions (non private scope only)

Name	Connector pane	Description	S.	R.	I.
setPosition		<p>!!! Attention: this VI can potentially uncalibrate your Mono !!! Only use after reading the description !!! Use moveToPosition to move the mono to a different wavelength !!!</p> <p>This command sets the wavelength value of the current grating position of the monochromator. This could potentially uncalibrate the monochromator and report an incorrect wavelength compared to the actual output wavelength.</p> <p>wavelength Float. Set the wavelength of the mono at the current position.</p>			
getFilterWheelPosition		<p>Returns the current filter wheel position.</p> <p>parameter description index Integer. Used to identify which mono to control. See mono_list command locationId Integer. Specifies the filter wheel location. 0 = Filter wheel 1 (Internal) 1 = Filter wheel 2 (External)</p>			
moveFilterWheel		<p>Move the filter wheel to a position.</p> <p>locationId Integer. Specifies which filter wheel to move. 0 = Filter wheel 1 (Internal) 1 = Filter wheel 2 (External) position Integer. Position to move the filter wheel.</p>			
getGratingPosition		<p>Returns the current grating turret position.</p> <p>Note: Prior to the initialization of the grating turret, this value may not reflect the actual position of the turret. To read the current position of the grating turret, please run mono_init prior to running this command.</p>			

Name	Connector pane	Description	S.	R.	I.
moveGrating		<p>Moves the grating turret to the specified position.</p> <p>Note: The turret sensor does not re-read the position each time it is moved, therefore the position may not be accurate prior to initialization. See note for mono_getGratingPosition.</p>			
getMirrorPosition		<p>Returns the position of the specified mirror.</p> <p>parameter description index Integer. Used to identify which mono to control. See mono_list command locationId Integer. Identifies which mirror to get the position from. 0 = Mirror 1 (Entrance) 1 = Mirror 2 (Exit)</p>			
moveMirror		<p>Moves the specified mirror to a position.</p> <p>locationId Integer. Identifies which mirror to move (zero-based). 0 = Mirror 1 (Entrance) 1 = Mirror 2 (Exit) position Integer. Position to move to. 0 = Axial 1 = Lateral</p>			
getPosition		<p>Returns the wavelength value, in nm, of the monochromator's current position.</p>			
moveToPosition		<p>This command starts the monochromator moving to the requested wavelength in nm. This is an asynchronous command. Use the mono_isBusy command to know when the move has completed.</p>			
getShutterStatus		<p>Returns the status of the currently selected shutter.</p> <p>Note: To view the status of the shutter solenoid the device must be configured for internal shutter mode.</p> <p>locationId Integer. Identifies the currently selected shutter. 0 = Shutter 1 (Front shutter) 1 = Shutter 2 (Side shutter) position Integer. Shutter position status. 0 = Closed 1 = Open</p>			

Name	Connector pane	Description	S.	R.	I.
shutterClose		<p>Deactivates the currently selected shutter solenoid.</p> <p>Note: The device must be configured for internal shutter mode. The shutter solenoid will not respond in External (Bypass) mode.</p>			
shutterOpen		<p>Activates the currently selected shutter solenoid.</p> <p>Note: The device must be configured for internal shutter mode. The shutter solenoid will not respond in External (Bypass) mode.</p>			
getSlitPositionInMM		<p>Returns the position of the specified slit in millimeters. The location id of each configured slit can be found under the ports section of the mono configuration. See mono_getConfig for additional information.</p> <p>For example:</p> <pre>"ports": [ { "locationId": 1, "slitType": 1 }, { "locationId": 2, "slitType": 1 }, { "locationId": 4, "slitType": 1 } ]</pre> <p>Note: The "locationId" parameter found in the mono configuration is 1-based. However, the mono_getSlitPositionInMM command uses a 0-based "locationId".</p>			

Name	Connector pane	Description	S.	R.	I.
moveSlitMM		<p>Moves the specified slit to the position in millimeters. The location id of each configured slit can be found under the ports section of the mono configuration. See mono_getConfig for additional information.</p> <p>For example:</p> <pre>"ports": [ { "locationId": 1, "slitType": 1 }, { "locationId": 2, "slitType": 1 }, { "locationId": 4, "slitType": 1 } ]</pre> <p>Note: The "locationId" parameter found in the mono configuration is 1-based. However, the mono_moveSlitMM command uses a 0-based "locationId".</p> <p>locationId Integer. Slit location (zero-based) position Float. Position in millimeters</p>			
DeviceClose		Closes communications with the monochromator indicated by the index.			
DeviceIsOpen		Returns true if selected monochromator is open.			
DeviceOpen		Opens communications with the monochromator indicated by the index command parameter.			
getConfig		<p>This command returns the monochromator configuration. Port Descriptions:</p> <p>locationId Integer. Used to identify the slit location. 1 = Front entrance (axial) 2 = Side entrance (lateral) 3 = Front exit (axial) 4 = Side exit (lateral) slitType Integer. Used to identify the slit size. 1 = 2mm slit 2 = 7mm slit</p>			

Name	Connector pane	Description	S.	R.	I.
Init		Starts the monochromator initialization process (homing...). This is a "long-running" asynchronous command. Use the mono_isBusy command to know when initialization has completed.  force Boolean. Force starts the initialization process.			
isBusy		Returns true if selected monochromator is busy.			
isInitialized		This command returns true when the mono is initialized. Otherwise it returns false.  Note: This command may also return false when the mono is busy with another command.			
Read DeviceTerm		Accessor VI for this class property.			
Write DeviceTerm		Accessor VI for this class property.			
Read DeviceType		Accessor VI for this class property.			
Write DeviceType		Accessor VI for this class property.			
Read Index		Accessor VI for this class property.			
Write Index		Accessor VI for this class property.			
Read productId		Accessor VI for this class property.			
Write productId		Accessor VI for this class property.			
Read serialNumber		Accessor VI for this class property.			
Write serialNumber		Accessor VI for this class property.			

Scope: → Protected | → Community

Reentrancy: → Preallocated reentrancy | → Shared reentrancy

Inlining:  → Inlined

## 3.5. CCD.lvclass

**Responsibility:** This class contains all functionality needed for Horiba's Multi-Channel-detectors.

**Version:** 1.0.0.3








Table 13. Functions (non private scope only)




Name	Connector pane	Description	S.	R.	I.
getSpeedAndGainSettings		Wrapper around the getSpeed and getGain VIs for convenience.			
parseConfigData		Helper VI to parse the configuration string to a LV cluster.			
parseImageData		Helper VI to parse the ccd data to a LV 2D image array.			
parseRangeMode		No description found (add content in vi description)			
parseSpectralData		Helper VI to parse the ccd data to a LV 1D image array and 1d Spectral information array.			
sortGetConfig		Helper VI to sort the config data.			
acquisitionAbort		Stops the current acquisition.			
acquisitionStart		Starts an acquisition that has been set up according to the previously defined acquisition parameters.  Note: To specify the acquisition parameters please see ccd_setROI and ccd_setXAxisConversionType. If there are no acquisition parameters specified at the time of acquisition it may result in no data being generated.			
calculateRangeModePositions		Starts an acquisition that has been set up according to the previously defined acquisition parameters.  Note: To specify the acquisition parameters please see ccd_setROI and ccd_setXAxisConversionType. If there are no acquisition parameters specified at the time of acquisition it may result in no data being generated.			
DeviceClose		Closes the connection to the CCD device.			
DeviceIsOpen		Checks if connection to the device is open or not.			
DeviceOpen		Opens the connection to the CCD device.			

Name	Connector pane	Description	S.	R.	I.
getAcqCount		<p>Gets the number of acquisition measurements to be perform sequentially by the hardware.</p> <p>Return Results:</p> <p>results description count Integer. The number of acquisition measurements to be performed.</p>			
getAcquisitionBusy		No description found (add content in vi description)			
getAcquisitionData		<p>The acquisition description string consists of the following information:</p> <p>acqIndex: Acquisition number roiIndex: Region of Interest number xOrigin: ROI's X Origin yOrigin: ROI's Y Origin xSize: ROI's X Size ySize: ROI's Y Size xBinning: ROI's X Bin yBinning: ROI's Y Bin Timestamp: This is a timestamp that relates to the time when the all the programmed acquisitions have completed. The data from all programmed acquisitions are retrieved from the CCD after all acquisitions have completed, therefore the same timestamp is used for all acquisitions. Command Parameters:</p> <p>Return Results:</p> <p>results description acquisition String. Acquisition data. Example command:</p> <p>Example response:</p> <pre>{ "command": "ccd_getAcquisitionData",   "errors": [], "id": 1234, "results": {     "acquisition": [ { "acqIndex": 1, "roi": [ {       "roiIndex": 1, "xBinning": 1, "xOrigin": 1,       "xSize": 8, "xyData": [ [ 885.6389770507812,         976 ], [w 885.2899780273438, 975 ], [         884.9409790039062, 979 ], [         884.593017578125, 976 ], }     ] }   ] }</pre>			
getAcquisitionReady		No description found (add content in vi description)			

Name	Connector pane	Description	S.	R.	I.
getChipSize		<p>Returns the chip sensor's pixel width and height size.</p> <p>Return Results:</p> <p>results description x Integer. Chip sensor's x size in pixels (width) y Integer. Chip sensor's y size in pixels (height)</p>			
getChipTemperature		<p>Returns the temperature of the chip sensor in degrees C.</p> <p>Return Results:</p> <p>temperature Float. Chip sensor temperature in degrees C.</p>			
getCleanCount		<p>Gets the number of cleans to be performed prior to measurement.</p> <p>results description count Integer. Number of cleans. mode Integer. Specifies how the cleans will be performed. 0 = Never 1 = First Only 2 = Between Only 3 = Each</p>			

Name	Connector pane	Description	S.	R.	I.
getConfig	 <p>CCD in [11] [3] CCD out error in (no error) [8] [2] config-JSON string Delay between send and rece... [6] [1] config data [0] error out</p>	<p>Returns the CCD device configuration.</p> <p>results description configuration JSON. CCD device configuration.</p> <p>example response:</p> <pre>{ "command": "ccd_getConfig", "errors": [],   "id": 1234, "results": { "configuration": {     "chipHSpacing": "140", "chipHeight": "70",     "chipName": "S10420",     "chipSerialNumber": "FAH23 098",     "chipVSpacing": "140", "chipWidth": "2048",     "deviceType": "HORIBA Scientific     Syncerity", "fitParameters": [ 0, 1, 0, 0, 0 ],     "gains": [ { "info": "Best Dynamic Range",       "token": 1 }, { "info": "High Sensitivity",       "token": 2 }, { "info": "High Light", "token": 0     } ], "hardwareAvgAvailable": false,     "lineScan": false, "parallelSpeeds": [ {       "info": "9.6 µSec", "token": 1 }, { "info": "4.9       µSec", "token": 2 }, { "info": "19 µSec",       "token": 0 } ], "productId": "13",     "serialNumber": "Camera SN: 5128",     "signals": [ { "events": [ { "name": "Ready       For Trigger", "token": 1, "types": [ { "name":       "TTL Active Low", "token": 1 }, { "name":       "TTL Active High", "token": 0 } ] }, { "name":       "Not Readout", "token": 2, "types": [ {       "name": "TTL Active Low", "token": 1 }, {       "name": "TTL Active High", "token": 0 } ] }, {       "name": "Shutter Open", "token": 3, "types":       [ { "name": "TTL Active Low", "token": 1 }, {       "name": "TTL Active High", "token": 0 } ] }, {       "name": "Start Experiment", "token": 0,       "types": [ { "name": "TTL Active Low",       "token": 1 }, { "name": "TTL Active High",       "token": 0 } ] } ], "name": "Signal Output",       "token": 0 } ], "speeds": [ { "info": "500 kHz       ", "token": 1 }, { "info": "500 kHz Ultra",       "token": 2 }, { "info": "500 kHz Wrap",       "token": 127 }, { "info": " 45 kHz ", "token": 0     } ], "supportedFeatures": {       "cf_3PositionSlit": false,       "cf_CMOSOffsetCorrection": false,       "cf_Cleaning": true, "cf_DSP": false,       "cf_DSPBin2X": false,       "cf_DelayAfterTrigger": false, "cf_Delays":</pre>			

Name	Connector pane	Description	S.	R.	I.
getDataSize		<p>Gets the number of pixels to be returned based on the current settings.</p> <p>results description size Integer. Byte data size for all ROIs and acquisitions.</p>			
getExposureTime		<p>Gets the exposure time (expressed in Timer Resolution units).</p> <p>Note: To check the current Timer Resolution value see <code>ccd_getTimerResolution</code>. Alternatively the Timer Resolution value can be set using <code>ccd_setTimerResolution</code>.</p> <p>Example: If Exposure Time is set to 50, and the Timer Resolution value is 1000, the CCD exposure time (integration time) = 50 milliseconds.</p> <p>If Exposure Time is set to 50, and the Timer Resolution value is 1, the CCD exposure time (integration time) = 50 microseconds.</p>			
getFitParams		<p>Gets the FIT parameters contained in the CCD configuration for the conversion of pixel to wavelength if done via the settings contained in the CCD.</p>			
getGain		<p>Gets the current gain token and the associated description information for the gain token. Gain tokens and their descriptions are part of the CCD configuration information. See <code>ccd_getConfig</code> command. For example:</p> <pre>"gains": [ { "info": "Best Dynamic Range", "token": 1 }, { "info": "High Sensitivity", "token": 2 }, { "info": "High Light", "token": 0 } ]</pre>			


Name	Connector pane	Description	S.	R.	I.
getParallelSpeed		<p>Gets the current parallel speed token and token description. Parallel speed tokens and their descriptions are contained in the CCD configuration information. See <code>ccd_getConfig</code> command.</p> <p>Note: The Parallel Speed value may also be referred to as the Vertical Shift Rate. These terms are interchangeable.</p> <p>For example:</p> <pre>"parallelSpeeds": [ { "info": "9.6 μSec", "token": 1 }, { "info": "4.9 μSec", "token": 2 }, { "info": "19 μSec", "token": 0 } ],</pre>			
getSignalOut		<p><code>ccd_getSignalOut</code> This command is used to get the current setting of the signal output. The address, event, and signalType parameters are used to define the signal based on the supported options of that particular CCD. The supported signal options are retrieved using the <code>ccd_getConfig</code> command, and begin with the "Signals" string contained in the configuration. For example:</p> <pre>"signals": [ { "events": [ { "name": "Shutter Open", "token": 3, "types": [ { "name": "TTL Active Low", "token": 1 }, { "name": "TTL Active High", "token": 0 } ] }, { "name": "Start Experiment", "token": 0, "types": [ { "name": "TTL Active Low", "token": 1 }, { "name": "TTL Active High", "token": 0 } ] } ], "name": "Signal Output", "token": 0 } ]</pre>			
getSpeed		<p><code>ccd_getSpeed</code> Gets the current speed token and the associated description information for the speed token. Speed tokens and their descriptions are part of the CCD configuration information. See <code>ccd_getConfig</code> command. For example:</p> <pre>"speeds": [ { "info": "500 kHz ", "token": 1 }, { "info": "500 kHz Ultra", "token": 2 }, { "info": "500 kHz Wrap", "token": 127 }, { "info": " 45 kHz ", "token": 0 } ]</pre>			

Name	Connector pane	Description	S.	R.	I.
getTimerResolution		<p>Gets the current timer resolution token.</p> <p>results description resolutionToken Integer. Timer resolution token. 0 - Timer resolution is set to 1000 microseconds 1 - Timer resolution is set to 1 microsecond</p>			
getTriggerIn		<p>This command is used to get the current setting of the input trigger. The address, event, and signalType parameters are used to define the input trigger based on the supported options of that particular CCD. The supported trigger options are retrieved using the ccd_getConfig command, and begin with the "Triggers" string contained in the configuration. For example:</p> <pre>"triggers": [ { "events": [ { "name": "Each - For Each Acq", "token": 1, "types": [ { "name": "TTL Rising Edge", "token": 1 }, { "name": "TTL Falling Edge", "token": 0 } ] }, { "name": "Once - Start All", "token": 0, "types": [ { "name": "TTL Rising Edge", "token": 1 }, { "name": "TTL Falling Edge", "token": 0 } ] } ], "name": "Trigger Input", "token": 0 } ]</pre>			
getXAxisConversionType		<p>Gets the X axis pixel conversion type to be used when retrieving the acquisition data with the ccd_getAcquisitionData command.</p> <p>results description type Integer. The X-axis pixel conversion type to be used. 0 = None (default) 1 = CCD FIT parameters contained in the CCD firmware 2 = Mono Wavelength parameters contained in the icl_settings.ini file</p>			
restart		<p>Performs a restart on the CCD.</p>			
setAcqCount		<p>Sets the number of acquisition measurements to be performed sequentially by the hardware. A count &gt; 1 is commonly referred to as "MultiAcq".</p>			






Name	Connector pane	Description	S.	R.	I.
setAcqFormat		<p>Sets the acquisition format and the number of ROIs (Regions of Interest) or areas. This command will remove all previously defined ROIs. After using this command, the ccd_setRoi command should be used to define each ROI.</p> <p>parameter description numberOfRois Integer. Number of ROIs (Regions of Interest / areas) format Integer. The acquisition format. 0 = Spectra 1 = Image 2 = Crop* 3 = Fast Kinetics* * Note: The Crop (2) and Fast Kinetics (3) acquisition formats are not supported by every CCD.</p>			
setCenterWavelength		<p>This command sets the center wavelength value and other parameters to be used in the pixel to wavelength conversion.</p> <p>Note: This command should be called before ccd_setXAxisConversionType and ccd_setAcquisitionStart and is only useful if the xAxisConversion type is set to Fitparams.</p>			
setCleanCount		<p>Sets the number of cleans to be performed according to the specified mode setting.</p> <p>parameter description index Integer. Used to identify which CCD to target. See ccd_list command count Integer. Number of cleans. mode Integer. Specifies how the cleans will be performed. 0 = Never 1 = First Only 2 = Between Only 3 = Each</p>			

Name	Connector pane	Description	S.	R.	I.
setExposureTime		<p>Sets the exposure time (expressed in Timer Resolution units).</p> <p>Note: To check the current Timer Resolution value see <code>ccd_getTimerResolution</code>. Alternatively the Timer Resolution value can be set using <code>ccd_setTimerResolution</code>.</p> <p>Example: If Exposure Time is set to 50, and the Timer Resolution value is 1000, the CCD exposure time (integration time) = 50 milliseconds.</p> <p>If Exposure Time is set to 50, and the Timer Resolution value is 1, the CCD exposure time (integration time) = 50 microseconds.</p>			
setGain		<p>Sets the CCD gain token. A list of supported gain tokens can be found in the CCD configuration. See <code>ccd_getConfig</code> command. For example:</p> <pre>"gains": [ { "info": "Best Dynamic Range", "token": 1 }, { "info": "High Sensitivity", "token": 2 }, { "info": "High Light", "token": 0 } ]</pre>			
setParallelSpeed		<p>Sets the CCD parallel speed token. A list of supported parallel speed tokens can be found in the CCD configuration. See <code>ccd_getConfig</code> command.</p> <p>Note: The Parallel Speed value may also be referred to as the Vertical Shift Rate. These terms are interchangeable.</p> <p>For example:</p> <pre>"parallelSpeeds": [ { "info": "9.6 μSec", "token": 1 }, { "info": "4.9 μSec", "token": 2 }, { "info": "19 μSec", "token": 0 } ],</pre>			

Name	Connector pane	Description	S.	R.	I.
setRoi		<p>Sets a single (roiIndex) ROI (Region of Interest) or area as defined by the X and Y origin, size, and bin parameters. The number of ROIs may be set using the ccd_setAcqFormat command. For Spectral acquisition format set yBin = ySize.</p> <p>Note: All values must fall within the x and y limits of the chip sensor, see ccd_getChipSize. If the ROI is not valid, the device will not be properly setup for acquisition.</p> <p>Command Parameters:</p> <p>parameter description index Integer. Used to identify which CCD to target. See ccd_list command roiIndex Integer. The region of interest's index (one-based) xOrigin Integer. The starting pixel in the x direction (zero-based) yOrigin Integer. The starting pixel in the y direction (zero-based) xSize Integer. The number of pixels in the x direction (one-based) ySize Integer. The number of pixels in the y direction (one-based) xBin Integer. The number of pixels to “bin” (x pixels summed to 1 value) yBin Integer. The number of pixels to “bin” (y pixels summed to 1 value)</p>			

Name	Connector pane	Description	S.	R.	I.
setSignalOut		<p>This command is used to enable or disable the signal output. When enabling the signal output, the address, event, and signalType parameters are used to define the signal based on the supported options of that particular CCD. The supported signal options are retrieved using the ccd_getConfig command, and begin with the "Signals" string contained in the configuration. For example:</p> <pre>"signals": [ { "events": [ { "name": "Shutter Open", "token": 3, "types": [ { "name": "TTL Active Low", "token": 1 }, { "name": "TTL Active High", "token": 0 } ] }, { "name": "Start Experiment", "token": 0, "types": [ { "name": "TTL Active Low", "token": 1 }, { "name": "TTL Active High", "token": 0 } ] }, { "name": "Signal Output", "token": 0 } ] }</pre> <p>parameter description index Integer. Used to identify which CCD to target. See ccd_list command enable Boolean. Enables or disables the signal. true = enable false = disable</p> <p>Note: When disabling the signal output, the address, event, and signalType parameters are ignored. address Integer. Token used to specify where the signal is located. (e.g. 0 = Signal Output)</p> <p>Note: Signal name and token can be found in the CCD config, see ccd_getConfig event Integer. Token used to specify when the signal event should occur. (e.g. 3 = Shutter Open)</p> <p>Note: Event name and token can be found in the CCD config, see ccd_getConfig signalType Integer. Token used to specify how the signal will cause the event. (e.g. 0 = TTL Active High)</p> <p>Note: Signal type and token can be found in the CCD config, see ccd_getConfig</p>			

Name	Connector pane	Description	S.	R.	I.
setSpeed		<p>Sets the CCD speed token. A list of supported speed tokens can be found in the CCD configuration. See ccd_getConfig command. For example:</p> <pre>"speeds": [ { "info": "500 kHz ", "token": 1 }, { "info": "500 kHz Ultra", "token": 2 }, { "info": "500 kHz Wrap", "token": 127 }, { "info": " 45 kHz ", "token": 0 } ]</pre>			
setTimerResolution		<p>Sets the current timer resolution token.</p> <p>resolutionToken Integer. Timer resolution token. 0 - Sets the timer resolution to 1000 microseconds 1 - Sets the timer resolution to 1 microsecond*</p>			
setTriggerIn		<p>This command is used to enable or disable the trigger input. When enabling the trigger input, the address, event, and signalType parameters are used to define the input trigger based on the supported options of that particular CCD. The supported trigger options are retrieved using the ccd_getConfig command, and begin with the "Triggers" string contained in the configuration. For example:</p> <pre>"triggers": [ { "events": [ { "name": "Each - For Each Acq", "token": 1, "types": [ { "name": "TTL Rising Edge", "token": 1 }, { "name": "TTL Falling Edge", "token": 0 } ] }, { "name": "Once - Start All", "token": 0, "types": [ { "name": "TTL Rising Edge", "token": 1 }, { "name": "TTL Falling Edge", "token": 0 } ] } ], "name": "Trigger Input", "token": 0 } ]</pre>			

Name	Connector pane	Description	S.	R.	I.
setXAxisConversionType		<p>Sets the X-axis pixel conversion type to be used when retrieving the acquisition data with the <code>ccd_getAcquisitionData</code> command.</p> <p>Note: To use the parameters contained in the <code>icl_settings.ini</code> file, the <code>ccd_setCenterWavelength</code> command must be called first.</p> <p>type Integer. The X-axis pixel conversion type to be used. 0 = None (default) 1 = CCD FIT parameters contained in the CCD firmware 2 = Mono Wavelength parameters contained in the <code>icl_settings.ini</code> file</p>			
Read DeviceTerm		After initialization of the device this property gives back the device term which is "ccd" for the device type CCD. This is used for prefixing all commands to the ICL.			
Write DeviceTerm		Sets the device term after init.			
Read DeviceType		Returns the device type of the device saved in the device firmware.			
Write DeviceType		Sets the device type after init			
Read Index		Returns the index of the device in the ICL layer.			
Write Index		Sets the device id after init			
Read productId		Returns the productID of the device saved in the device firmware.			
Write productId		Sets the device productID after init			
Read serialNumber		Returns the serial of the device saved in the device firmware.			
Write serialNumber		Sets the device serial after init			
Setup Example		No description found (add content in vi description)			

Name	Connector pane	Description	S.	R.	I.
Teardown Example	Module Was Already Running? [11]  [0] error out	No description found (add content in vi description)			
Example - Open Close CCD	CCD in [11]  [3] CCD out [2] device is open? [0] error out	No description found (add content in vi description)			
Example - Get Config	CCD in [11]  [3] CCD out [2] config data [0] error out	No description found (add content in vi description)			
CCD Example	exposure time [ms] [0]  [4] ExposureTime [ms] [6] device is open? [15] error out [14] chip size, y [13] chip size, x	Example for a CCD spectrum acquisition.			

Scope: → Protected | → Community

Reentrancy: → Preallocated reentrancy | → Shared reentrancy

Inlining: → Inlined

## 3.6. SCD.lvclass

**Responsibility:** This class contains all functionality needed for Horiba's Single-Channel-detectors.

**Version:** 1.0.0.1




Table 14. Functions (non private scope only)

Name	Connector pane	Description	S.	R.	I.
DeviceClose	SCD in [11]  [3] SCD out [2] results-JSON string [1] ID [0] error out [4] device is open?	Closes communications with the CCD indicated by the index.			
DeviceOpen	SCD in [11]  [3] SCD out [2] device is open? [1] ID [0] error out	Opens communications with the SpectrAcq3 indicated by the index command parameter.			
DeviceIsOpen	SCD in [11]  [3] SCD out [2] device is open? [1] ID [0] error out	Returns true if selected SpectrAcq3 is open.			
Read DeviceTerm	SCD in [11]  [3] SCD out [2] deviceTerm [0] error out	Accessor VI for this class property.			

Name	Connector pane	Description	S.	R.	I.
Write DeviceTerm		Accessor VI for this class property.			
Read DeviceType		Accessor VI for this class property.			
Write DeviceType		Accessor VI for this class property.			
Read Index		Accessor VI for this class property.			
Write Index		Accessor VI for this class property.			
Read productId		Accessor VI for this class property.			
Write productId		Accessor VI for this class property.			
Read serialNumber		Accessor VI for this class property.			
Write serialNumber		Accessor VI for this class property.			

Scope:  → Protected |  → Community

Reentrancy:  → Preallocated reentrancy |  → Shared reentrancy

Inlining:  → Inlined



# Chapter 4. Custom errors

**TIP** Custom errors are added via vi named `*--error.vi`.

Table 15. Custom errors

Name	Code	Description	Owned by
Module Not Running	0		<a href="#">DeviceManager.lvlib</a>
Module Not Stopped	0		<a href="#">DeviceManager.lvlib</a>
Module Not Synced	0		<a href="#">DeviceManager.lvlib</a>
Request and Wait for Reply Timeout	0		<a href="#">DeviceManager.lvlib</a>

# Chapter 5. Legal Information

## 5.1. Document creation

This document has been generated using the following tools.

### 5.1.1. Antidoc

Project website: [Antidoc](#)

Maintainer website: [Wovalab](#)

BSD 3-Clause License

Copyright © 2019, Wovalab, All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
- Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
- Neither the name of the copyright holder nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT HOLDER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

### 5.1.2. AsciiDoc for LabVIEW™

Project website: [AsciiDoc toolkit](#)

Maintainer website: [Wovalab](#)

BSD 3-Clause License

Copyright © 2019, Wovalab, All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
- Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
- Neither the name of the copyright holder nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT HOLDER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

### 5.1.3. Graph Builder

Project website: [Graph Builder](#)

BSD 3-Clause License

Copyright © 2020, Cyril GAMBINI All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
- Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
- Neither the name of the copyright holder nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT HOLDER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES

(INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

#### 5.1.4. **classy Diagram Viewer**

Project website: [classy Diagram Viewer](#)

BSD 3-Clause License

Copyright © 2021, Tatiana Boyé All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
- Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
- Neither the name of the copyright holder nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT HOLDER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

## 5.2. Product used in the project

The documented project has been developed with the following products.

### 5.2.1. DQMH®

Copyright © 2021 DQMH® Consortium, LLC. All Rights Reserved.

Find more details on [DQMH® Consortium](#) website